



Malaria P.f/ PAN Rapid Card

Cat.No	Product	Type (pack size)	Specimen	Shelf life	Storage
1/012	Malaria Pf/PAN	Card (301/KIT)	Whole Blood	24M	1-40°C

General Information

Malaria Pf/ PAN, Pf/P.v, Pf

Pf/PAN Sensitivity: 99.6% (Pf) 95.9% (Non-Pf)

Pf/P.v Sensitivity: 99.6% (Pf) 98% (P.v)

Specificity: 99.3%

Specimen: Whole blood 5µl using loop

Sample diluent: 4 drops (120µl)

Reading Time: 20 min