



Malaria P.f/ PAN Rapid Card

Cat.No	Product	Type (pack size)	Specimen	Shelf life	Storage
17012	Malaria P.f/PAN	Card (30T/Kit)	Whole Blood	24M	1~40°C

General Information

Malaria P.f/ PAN, P.f/P.v, P.f

P.f/PAN Sensitivity: 99.6% (P.f) 95.9% (Non-P.f)

P.f/P.v Sensitivity: 99.6% (P.f) 98% (P.v)

Specificity: 99.3%

Specimen: Whole blood 5µl using loop

Sample diluent: 4 drops (120µl)

Reading Time: 20 min